



Linking Users With Hearing Experts

A Cloud-Based Solution For Remote Fitting

Overview

Jacoti earCloud® is the web service that allows users to store their audio profiles in the cloud, share them across their devices and connect to the Hearing Expert of their choice for remote fitting*.


Jacoti earCloud® is included in the Jacoti Inside solution integrated in Qualcomm QCC5100 SoCs and available on Cadence's Tensilica HiFi digital signal processors as cloud and backend solution.

* Remote fitting is the process of performing fitting without the physical presence of the hearing expert.

- ▶ Secure web service
- ▶ Full integration with Jacoti ListenApp®, Jacoti Hearing Center and Jacoti Hearing Center Pro
- ▶ Cloud and backend solution for Jacoti Inside
- ▶ Hearing profile synching across devices
- ▶ Hearing expert linking

Regulatory and clinical information

- ▶ Scalable and Secure medical storage cloud
- ▶ Possibility to remove all user details and/or medical data as well any subscriptions to newsletter



Software components

Linking with hearing experts as a user	The user can choose and invite their Hearing Expert
Remote fitting as an expert	The expert can link to the user by knowing their email address and perform remote fitting on their hearing loss
User audiogram manual input	The expert can manually edit the audiogram of the user
User clinical history	The expert can have access to all the hearing test results of the user thus having a clear picture of how he or she evolves.

Jacoti earCloud® is accessible from a web browser at earcloud.net
Microsoft Edge 90 and above, recent Chrome and Mozilla Firefox versions among others.

About Jacoti

Jacoti BV | Hearing Technologies is a science-based company that develops hearing enhancement solutions embeddable in consumer devices. Its flagship product, Jacoti Inside, optimizes audio to each individual hearing requirement from consumer technologies to fully-fledged medical devices.

Contact us: www.jacoti.com | info@jacoti.com